

# MUTARS®

---



implantcast



## Distal Femur M-O-M

### surgical technique

Femorotibial M-O-M coupling





# Distal Femur M-O-M

## surgical technique

Femorotibial M-O-M coupling

MUTARS® was developed in co-operation with Prof. Dr. W. Winkelmann (former director) and Prof. Dr. G. Gosheger (director), Clinic and Polyclinic for General Orthopedics and Tumororthopedics at the University Hospital of Münster, Germany. MUTARS® has been in successful clinical use since 1992.

## Table of Contents

Table of Contents .....	1
Information about the Silver and TiN coating .....	2
System overview.....	3
Assembling options/reconstruction length.....	4
Surgical technique .....	5
Implants .....	23
Instruments .....	27

**Nota Bene:** The described surgical technique is the suggested treatment for the uncomplicated procedure. In the final analysis the preferred treatment is that which addresses the needs of the individual patient.

**Copyright note:** MUTARS®, implavit® and implatan® are registered trademarks of implantcast GmbH. The use and copy of the content of this brochure are only allowed with prior permit given by the implantcast GmbH.



### **The Silver coating**

Early and late infections represent the most severe complications of tumour arthroplastic treatments. Although local and systemic antibiotic treatments are considered, the scientific literature reports of infection rates from 5 to 35 percent. Reasons for these high rates are, for example, the long surgery time, the large incisions and the immunosuppression due to chemo therapy and radio therapy as well as the increasing resistance of the bacteria against antibiotic drugs.

The anti-infective effect of silver ions has been known for centuries i.e. the disinfection of potable water is based on this principle. This special property of silver is used for the silver coated components of MUTARS<sup>®</sup> to build an intelligent protection against bacteria.

Until now only non-articulating surfaces and surfaces without direct bony contact are coated with silver.

In the catalogue information of this brochure you can find the supplement \*S indicating which MUTARS<sup>®</sup> components are available in a silver coated version. The eight digit REF number receives an addition after the last digit (e.g. 5220-0020S).

### **It is not permitted to flush the wound with antiseptics that contain iodine or heavy metals (such as Betaisodona<sup>®</sup>)**

Iodine and Silver form insoluble salt complexes not only with the silver ions that are released post-operatively but also with the silver layer of the implant that will be covered with an insoluble silver-iodine (AgI) film. This will destroy the anti-adhesive protective layer irreversibly. Iodine or heavy metal based antiseptics may not be used at any time. Alternatively solutions containing H<sub>2</sub>O<sub>2</sub> – (like Lavasept<sup>®</sup>, Prontosan<sup>®</sup> or similar) can be used.

The silver coating can be destroyed in its function by two factors: large amounts of albumin from seroma or hematoma can bind larger amounts of silver (1 mol Albumin inactivates 3 moles Silver ions). This should be minimized by using an attachment tube. In the instance that an infection is known pre-operatively, antibiotics like Vancomycin can be mixed with the bone cement. The intramedullary stems are not silver coated and cemented components are preferred in case of a septic revision.

### **The TiN coating for allergy prophylaxis**

As the metallic components of total knee replacements, the articulating metallic parts of the MUTARS<sup>®</sup> system are made of casted CoCrMo alloy. In the late 70's and 80's of the last century, some of the Cobalt Chromium implants had a small Nickel content to add strength to the implant. Nickel is the primary cause for metal sensitivity, although some patients have shown to be hypersensitive to other metals such as Cobalt and Chromium. The use of titanium components can't solve this problem, because the wear of the articulating polyethylene inlays will increase and so the survival time of the prosthesis is reduced. Since the end of the 1990's TiN (Titanium Nitride coating) has been successfully applied to protect the body against metal ions that could cause allergic reactions.

The metal ion release of TiN coated or TiNbN coated implants is reduced down to 10%.<sup>1</sup>

In order to prevent allergic reactions, certain parts of the prosthesis may be supplied with a ceramic coating (TiN). Since almost all components of the tumor system consist of titanium alloy, this only concerns those components, which are made of a cast CoCr alloy (CoCrMo). The REF-numbers of the TiN coated implants have the suffix N after the last digit (e.g. 5720-0005N).


\*S: For anti-infective treatment, silver coated implants are available.

\*N: For anti-allergic treatment, TiN coated implants are available.

\*SN: Implants are coated with silver and TiN.

<sup>1</sup> Metal Ion Release from Non-Coated and Ceramic Coated Femoral Knee Components: Boil test 240h in NaCl-solution nach FMZ PhysWerk VA 97350, University Würzburg (D) (On File)

# System Overview





# MUTARS<sup>®</sup> Distal Femur

## distal femur replacement assembling options (length in mm)

reconstruction	distal femur	components		
		connecting part 100 mm	extension piece	bar screw
100 mm	90*			25
120 mm	110	-	-	45
140 mm	90*		40	65
160 mm	110	-	40	85
180 mm	110	-	60	105
200 mm	110	-	80	125
220 mm	110	100	-	45 + 25
240 mm	110	-	80 + 40	165
260 mm	110	100	40	65 + 45
280 mm	110	100	60	85 + 45
300 mm	110	100	80	105 + 45
320 mm	110	100	60 + 40	125 + 45

\*A distal femur 90 mm is available on special request  
(reconstruction length 100 mm)

**Note:** Please notice that the amount of implants and instruments send with an individual shipment may differ from the information in the catalogue information of this brochure. Please make sure, during the preoperatively planning, that all necessary implants and instruments are available for the surgery.




figure 1a and 1b




figure 2a and 2b

## Tumor resection

Resect the tumour and measure the length of the explant.

The minimum bone resection should be 120mm (or 100 mm if the special distal femur 90 mm is used, available on demand). Remove the menisci.

## Tibial bone preparation

Open the tibial medullary cavity with the universal drill  $\varnothing$  6 mm (fig. 1a and 1b). The drilling should be orientated to open the center of the medullary cavity (eminentia intercondylaris: ventral 1/3, dorsal 2/3).

Enlarge the opening of the medullary cavity with rigid drills (fig. 2a and 2b).

To choose the correct reamer size for the use of a **cementless tibial stem** consult table 1, for the use of a **cemented tibial stem** consult table 2.

**Table 1: cementless implantation**

Tibial stem 12 mm	→ drill 11 mm
Tibial stem 14 mm	→ drill 13 mm
Tibial stem 16 mm	→ drill 15 mm
Tibial stem 18 mm	→ drill 17 mm

**Table 2: cemented implantation**

Tibial stem 11 mm	→ drill 13 mm
Tibial stem 13 mm	→ drill 15 mm
Tibial stem 15 mm	→ drill 17 mm
Tibial stem 17 mm	→ drill 17 mm

To ascertain adequate depth is met, the drills have depth marks (120 mm for 120 mm stems, 160 mm for 160 mm stems and 200 mm for 200 mm stems) corresponding with the tibial stem length (fig. 2a and 2b). The last drill used is left in the tibial canal.

The tibia resection block 0° is attached to the intramedullary tibial alignment guide and the cutting block is placed over the tibial drill left in the intramedullary canal (fig. 3a).

Adjust the rotational alignment and lock the alignment guide by impacting the two spikes into the tibial surface and lock all quick connectors (fig. 3b).

Slide the tibial stylus into the upper slot of the resection block to adjust the resection height. Make sure that the marking SLOTTED<sub>1</sub> is directed to the bone, when a slotted cut is planned (fig. 4). If a nonslotted cut should be performed the NONSLOTTED marking on the stylus should point to the bone.

For the primary bone cut, make sure that the stylus is adjusted to the 15 mm mark<sub>2</sub> and 15mm of bone will be removed from the tibia (fig. 4).

In revision cases normally a minimum bone cut is recommended and the stylus should be adjusted to the 2mm height. When the correct resection height is determined, please lock the quick connector at the resection block.

Please insert the fixation pins in the marked level to fix the block to the bone. Remove the tibial resection stylus. If necessary please use the 3,2mm drill to predrill the holes (fig. 5).



and 3b




figure 4



figure 5

## MUTARS<sup>®</sup> Distal Femur M-O-M

---

Double check the resection angle and height by using the resection check (fig. 6a).



**figure 6a**



**figure 6b**

Use the ACS<sup>®</sup> saw blade to resect the bone. Prevent damaging of the intramedullary drill. If necessary please remove the drill before resectioning. For additional stability a pin can inserted in the oblique hole (fig. 6b).

Please check the quality of the cut. Make sure that the cut is totally flat and remove the resection block.



**figure 7**

The resected tibia is checked and the reamer guide with the tibial centering guide is slide over the tibial reamer in place (fig. 7). The mark MEDIAL should be placed correctly to the medial side.



The right mediolateral alignment should be established and the tibial reamer guide is fixed with two pins (fig. 8a). The tibia reamer guide and the intramedullary tibial reamer are now removed.



figure 8a      figure 8b

### Use of tibial spacer

The joint line can be restored using tibial spacers or bone grafts. If necessary, additional bone should be resected to accommodate the trial tibial spacer. The trial tibial spacer is clicked under the tibial reamer guide. (fig. 8b). The height of the spacer should correspond with the one fixed at the preoperative assessment.

Combine the tibial reamer and the T-handle and ream carefully until the reamer is stopped by the chimney of the reamer (fig. 9a and 9b). It is strongly recommended **not** to use power tools for the reaming.



figure 9a      figure 9b

# MUTARS<sup>®</sup> Distal Femur M-O-M

The tibial fin punch is used to continue the tibial preparation. The punch should be punched down until it is stopped by the tibial reamer guide (fig. 10a and 10b).



figure 10a figure 10b

In case of sclerotic bone the tibial drill can be used. A drill sleeve is placed inside the tibia reamer guide to accommodate this drill. The drill sleeve is placed medially and the canal is drilled. After turning the sleeve 180°, the lateral side is to be drilled.

Remove all instruments.



figure 11a figure 11b

## Femoral bone preparation

Prepare the femoral medullary cavity with the MUTARS<sup>®</sup> medullary cavity reamer (fig. 12).



figure 12

## Cementless fixation

Ream the femoral medullary cavity preferably up to a depth of 130 mm with a flexible reamer 1,5 mm smaller than the preoperatively chosen femoral stem (fig. 13).

## Cemented fixation

Ream the femoral medullary cavity preferably up to a depth of 130 mm with a flexible reamer 2 mm larger than the preoperatively chosen femoral stem (fig. 13).

## Remark

In case no flexible reamers are in the hospital's stock flexible reamers can be provided on special demand.



figure 13



figure 14




figure 15

## Rasping of the femoral cavity

Assemble the femur rasp of the appropriate size (see tables below), the sleeve and the slide hammer. Lock the rasp on the slide hammer by using the engineers' wrench.

### Remark

The use of a femoral rasp for a **cemented stem** is optional. Generally you can proceed with the trial reduction (see page 11).

Mark the anterior aspect of the femoral bone to meet the correct antecurvatur of the femur (fig. 15).

### Use of cementless stems

Use the femur rasp (fig. 15), of the same size as the preoperatively chosen femoral stem (table 3).

Stem size	Rasp size
12mm	12mm
13mm	13mm
14mm	14mm
15mm	15mm
16mm	16mm
17mm	17mm
18mm	18mm

table 3

Rasp the medullary cavity with the chosen femoral rasp (fig. 15). A carefully use of the slide hammer is recommended.



## Optional technique for the use of cemented stems

If you want to prepare for a cemented stem with the femoral rasp, please use the rasp which is 2 mm larger than the preoperatively chosen cemented femoral stem (fig. 16).

That will provide a cement mantle of 1 mm thickness (table 4). Use the 18 mm rasp to prepare for the 17 mm stem.

Stem size	Rasp size
11mm	13mm
13mm	15mm
15mm	17mm
17mm	18mm

table 4

### Remark

It is recommended to clean the rasp of bone chips during the rasping.

Leave the femoral rasp in the bone for the trialing.




figure 16




figure 17

## Trial reduction

Mount the MUTARS<sup>®</sup> Distal Femur and the possibly needed extension pieces (possible enlargement from 20 to 260 mm; see table page 2) to the top of the rasp (fig. 17).

## Remark

For the **cemented procedure** bone rasps are usually not available. Please insert the cemented stem (without cement) or the trial stem for trialing purposes.

At that stage the use of a screw is optional, because the teeth mechanism gives the assembly a reasonable stability (fig. 18a and 18b).



figure 18a and 18b

## Use of trial implants

To check all resections performed, the tibial and femoral trial implants are used.

Screw the trial stem under the tibial trial of the selected size (fig. 19).

The stem is medialized and care should be taken to place the trial stem into the correct medio-lateral position. If necessary a trial spacer can be clicked under the trial tibial implant (fig. 19).



figure 19

The tibial trial and stem can be inserted using the tibial impactor (fig. 20a and 20b).



figure 20a figure 20b

The corresponding trial inlay is then placed on the tibial trial implant (fig. 21a), using the PE-inlay setting instrument (fig. 21b).




figure 21a figure 21b

Perform a trial reduction to assure that the correct femoral rotation is achieved and the joint line is restored in the correct height (fig. 22).

### Remark

Please notice that a coupling of the joint components is not possible at this stage when using the trial components. At a later stage the joint stability can be checked using the final implant components and the locking mechanism.



figure 22



figure 23a



figure 23b

Remove the trial inlay. Screw the slide hammer into the tapered hole of the trial plateau and remove the trial components (fig. 23a and 23b).

## Tibial component assembly



figure 24

Attach the selected tibial stem onto the cone of the tibia component and connect the two parts with the screws provided. An torque wrench 3.5 mm hex screw driver (fig. 24) should be used. The same way any tibial spacers should be added (fig. 24).



## MUTARS® Distal Femur M-O-M

Impact the tibial components with the tibial impactor (fig. 25a and 25b).



figure 25a figure 25b

After cement hardening, insert the PE-Inlay in the tibial joint. Insert the inlay from behind, move it forward towards the anterior locking rim and push it down at the posterior part until it is locked securely (fig. 26a). Consider to use the impactor for PE-Inlay (fig. 26b).

Although trial inserts are available, it is recommended to insert the final PE-inlay at that time in order to reduce the surgery time.



figure 26a



figure 26b

## Implantation of the femoral stem

Impact the MUTARS<sup>®</sup> femoral stem (fig. 27).

Insert the stem of the same size as the rasp if a **cementless stem** is used.

To prevent fractures of the cortical bone it is helpful to fix a bone forceps around the femoral bone during impaction.

If a cemented implantation is planned, insert the cement and use the **cemented stem** which is 2 mm smaller than the previously used reamer or rasp.

Remove all instruments during the cement hardening to prevent bending moments.




figure 27



### Mounting of the distal femoral implant components

Combine the Distal Femur and possibly needed extension pieces with the femoral stem. Make sure that the correct rotation of the distal femur is achieved. Insert the bar screw of the correct length (see table on page 2) (fig. 28a and 28b).

Lock the screw with the swing wrench while countering the assembly with the engineers' wrench (fig. 28b).

Insert the safety screw and lock it in the same way (fig. 28c).




figure 28b



figure 28a



figure 28c

## MUTARS<sup>®</sup> Distal Femur M-O-M

Assemble the MUTARS<sup>®</sup> locking mechanism and the special MUTARS<sup>®</sup> instrument for locking mechanism. Therefore turn the attachment part of the lock by 100 degrees until it rests in the sleeve of the locking instrument (fig. 29).




figure 29

Insert the lock into the intracondylar notch of the femoral joint (fig. 30).



figure 30

Use the socket wrench to turn the locking instrument and the lock clockwise by 180 degrees (fig. 31).



figure 31



## MUTARS® Distal Femur M-O-M

The lock is correctly positioned when the attachment part falls out of the sleeve of the locking instrument (fig. 32a). Remove the locking instrument.

The instrument to insert the mechanism into the tibia component is now placed in the top hole<sub>1</sub> of the coupling and the mechanism is guided into the hole of the tibial plateau (fig. 32b).

The coupling mechanism should be fully engaged and placed in the correct rotational position. The screw hole<sub>2</sub> (fig. 32a) should be placed forward-turned to enable locking (fig. 32b).

The positioner is inserted into the screw hole of the short stem of the coupling mechanism (fig. 33a and 33b).

The coupling mechanism is held in place with the setting instrument while removing the positioner. Then the locking bolt is fixed with the torque wrench 3.5 mm hex screw driver into the tibial component (fig. 34a).

Be sure that the locking bolt is fully engaged into the tibia component (fig. 34b).




figure 32b




figure 33b

figure 34a



figure 34b

## MUTARS<sup>®</sup> Distal Femur M-O-M

---



**figure 35a**

A hexagonal torque wrench 3.5 mm hex screw driver is used to screw in the Multi-lock security screw and to tighten the locking bolt (fig. 35a and 35b).



**figure 35b**



**figure 36a**

The implantation of the implant is now concluded. Stability and range of motion should be performed in flexion (fig. 36a) and extension (fig. 36b).



**figure 36b**

## Removal of an implant

In case a tibia component should be removed the Multilock security screw and the locking bolt should be removed from ventrally using the torque wrench 3,5mm hex screw driver.

The locking instrument is then used to remove the locking mechanism from the femoral component.

The femoral component can now be removed using the slide hammer and the special extractor (fig. 37a and 37b).

The tibial extractor is now attached to the slide hammer and placed into the screw hole of the tibial component (fig. 37a).

The attachment is secured using the rod with the small chain (fig. 37b).

The tibial component is now removed using the slide hammer (fig. 37c).



figure 37a



figure 37b



figure 37c

# MUTARS® Distal Femur M-O-M

## IMPLANTS

**\*S:** For anti-infective treatment, silver coated implants are available.

**\*N:** For anti-allergic treatment, TiN coated implants are available.

**\*SN:** Implants with Silver and TiN coating!



## MUTARS® Distal Femur M-O-M, incl. safety screw \*S \*N \*SN

*mat.: implavit®; CoCrMo-casting alloy according to DIN ISO 5832/4*

5720-0045	110 mm	left
5720-0040	110 mm	right
5720-0047	90 mm	left
5720-0042	90 mm	right



## MUTARS® extension piece \*S

*mat.: implatan®; TiAl<sub>6</sub>V<sub>4</sub> according to DIN ISO 5832/3*

5772-2504	40 mm
5772-2506	60 mm
5772-2508	80 mm
5772-2510	100 mm



## MUTARS® connecting part \*S

*mat.: implatan®; TiAl<sub>6</sub>V<sub>4</sub> according to DIN ISO 5832/3*

5730-0100	100 mm
-----------	--------



## MUTARS® attachment tube

*mat.: polyethylenterephthalat*

5900-0300	35 mm
5900-0310	55 mm



## IMPLANTS

### MUTARS® screw

mat.: *implatan*®; TiAl<sub>6</sub>V<sub>4</sub> according to DIN ISO 5832/3

5792-1002	M10x 25 mm
5792-1004	M10x 45 mm
5792-1006	M10x 65 mm
5792-1008	M10x 85 mm
5792-1010	M10x105 mm
5792-1012	M10x125 mm
5792-1014	M10x145 mm
5792-1016	M10x165 mm
5792-1018	M10x185 mm
5792-1020	M10x205 mm
5792-1022	M10x225 mm

### MUTARS® femoral stem cemented \*N

mat.: *implavit*®; CoCrMo-casting alloy according to DIN ISO 5832/3

5760-0011	11 mm
5760-0013	13 mm
5760-0015	15 mm
5760-0017	17 mm

### MUTARS® femoral stem cementless

mat.: *implatan*®; TiAl<sub>6</sub>V<sub>4</sub> according to DIN ISO 5832/3 with HA-coating

5760-0012	12 mm
5760-0113	13 mm
5760-0014	14 mm
5760-0115	15 mm
5760-0016	16 mm
5760-0117	17 mm
5760-0018	18 mm

### MUTARS® tibial plateau M-O-M \*N cementless, incl. screw for locking mechanism and safety screw

mat.: *implavit*®; CoCrMo-casting alloy acc. to DIN ISO 5832/4 Screw *implatan*®; TiAl<sub>6</sub>V<sub>4</sub> acc. to DIN ISO 5832/3 with TiN coating

5751-0203	xsmall
5751-0200	small
5751-0205	standard
5751-0210	large





## IMPLANTS

### MUTARS® tibial plateau M-O-M \*N cemented, incl. screw for locking mechanism and safety screw

mat.: *implavit*®; CoCrMo-casting alloy acc. to DIN ISO 5832/4 Screw *implatan*®; TiAl<sub>6</sub>V<sub>4</sub> acc. to DIN ISO 5832/3 with TiN coating

5751-0303	xsmall
5751-0300	small *S
5751-0305	standard *S
5751-0310	large *S

### screw for locking mechanism

mat.: *implatan*®; TiAl<sub>6</sub>V<sub>4</sub> acc. to DIN ISO 5832/3 with TiN coating

5720-1201

### MUTARS® screws for tibial plateau M-O-M (2 pcs)

mat.: *implatan*®; TiAl<sub>6</sub>V<sub>4</sub> according to DIN ISO 5832/3

5720-1205

### MUTARS® tibial stem, cementless

mat.: *implatan*®; TiAl<sub>6</sub>V<sub>4</sub> according to DIN ISO 5832/3

5756-1212	12 x 120 mm
5756-1214	14 x 120 mm
5756-1216	16 x 120 mm
5756-1218	18 x 120 mm
5756-1612	12 x 160 mm
5756-1614	14 x 160 mm
5756-1616	16 x 160 mm
5756-1618	18 x 160 mm
5756-2012	12 x 200 mm
5756-2014	14 x 200 mm
5756-2016	16 x 200 mm
5756-2018	18 x 200 mm

### MUTARS® tibial stem, cemented \*N

mat.: *implavit*®; CoCrMo-casting alloy according to DIN ISO 5832/4

5755-1211	11 x 120 mm
5755-1213	13 x 120 mm
5755-1215	15 x 120 mm
5755-1611	11 x 160 mm
5755-1613	13 x 160 mm
5755-1615	15 x 160 mm
5755-2011	11 x 200 mm
5755-2013	13 x 200 mm
5755-2015	15 x 200 mm



## IMPLANTS

### MUTARS® PE-inlay

mat.: UHMWPE according to

DIN ISO 5834/2

5721-0013 xsmall

5721-0002 small

5721-0001 standard

5721-0006 large



### MUTARS® locking mechanism

#### M-O-M

mat.: implavit®; CoCrMo-casting alloy according to

DIN ISO 5832/4

5720-1200



### MUTARS® patellar component, cemented

Mat.: UHMW-PE acc. to DIN ISO 5834/2

5720-1000



### Intramedullary plug

Mat.: UHMW-PE acc. to DIN ISO 5834/2

0299-4000 small

0299-4010 large



### MUTARS® tibial spacer, unicondylar \*S

Mat.: implatan®; TiAl<sub>6</sub>V<sub>4</sub> nach DIN ISO 5832/3

5810-0500 5 mm rl/lm

5810-1000 10 mm rl/lm.

5810-1500 15 mm rl/lm.

5810-2000 20 mm rl/lm.

5805-0500 5 mm ll/rm

5805-1000 10 mm ll/rm.

5805-1500 15 mm ll/rm

5805-2000 20 mm ll/rm.



### MUTARS® tibial spacer small \*S

Mat.: implatan®; TiAl<sub>6</sub>V<sub>4</sub> acc. to DIN ISO 5832/3

5800-2500 25 mm small right/left

5800-3505 35 mm small left

5800-5005 50 mm small left

5800-3500 35 mm small right

5800-5000 50 mm small right



### MUTARS® screw for tibial spacer M-O-M

Mat.: implatan®; TiAl<sub>6</sub>V<sub>4</sub>

according to DIN ISO 5832/3

5720-1203 for 5mm spacers

5720-1204 for 10-50 mm spacers







## MUTARS® Distal Femur M-O-M

---


### INSTRUMENTS




**MUTARS® basic instrument tray**  
7999-5712



**MUTARS® modular tibia instrument container I**  
7999-5733



**MUTARS® modular tibia instrument container II**  
7999-5738



**MUTARS® tibia modular trial container**  
7999-5736




## MUTARS<sup>®</sup> Distal Femur M-O-M


---

### INSTRUMENTS

**MUTARS<sup>®</sup> trial component tray**  
7999-7701



**MUTARS<sup>®</sup> rigid drills**  
7999-5735



**MUTARS<sup>®</sup> distal femoral M-O-M trial component tray**  
7999-7733 left  
7999-7734 right



# MUTARS® Distal Femur M-O-M

## INSTRUMENTS

### Content MUTARS® basic container



**MUTARS® universal impactor**  
7210-0000



**MUTARS® impact and extract sleeve**  
7230-0000



**MUTARS® socket wrench**  
7420-0000



**MUTARS® swing wrench**  
7411-0000



**MUTARS® engineers wrench SW 24**  
7490-0000



**MUTARS® slide hammer**  
7220-0001



**MUTARS® rasp for femoral stem**  
7760-0112                    12 mm  
7760-0113                    13 mm  
7760-0114                    14 mm  
7760-0115                    15 mm  
7760-0116                    16 mm  
7760-0117                    17 mm  
7760-0118                    18 mm



**handle for intramedullary plug**  
7512-4001



**MUTARS® medullary cavity reamer**  
7760 0501



## INSTRUMENTS

### Content MUTARS Tibia Container 1

**MUTARS<sup>®</sup> tibial impactor**  
7800-0008

**MUTARS<sup>®</sup> spacer block**  
7755-0010

**MUTARS<sup>®</sup> spacer block rotation**  
7755-0023

**hexagonal screw driver ¼" chuck, 3.5 mm**  
7512-0009

**torque limiter ¼" chuck 7Nm**  
7512-0007

**tibia cutting block revision 0°**  
7755-0054

**I / M tibial alignment guide**  
7755-0024

**MUTARS<sup>®</sup> tibial reamer guide**  
77550025 ± 2,5 ap  
7755-0039 ± 2,5 ap x-small

**fixation pin: 77 mm, D: 3,2 mm**  
4223-0029

**Universal drill 6 mm**  
7630-0106





# MUTARS® Distal Femur M-O-M

---

## INSTRUMENTS



**tibial resection stylus 15 mm**  
7700-0415



**MUTARS® tibial centralizer sleeve  
5mm**  
7755-0008



**MUTARS® sleeve for tibial preparator**  
7755-0022



**MUTARS® patella drill**  
7351-0000



**MUTARS® trial inlay**  
7721-0013      xsmall  
7721-0001      standard  
7721-0002      small  
7721-0006      large



**MUTARS® impactor for PE-inlay**  
7210-0001



## MUTARS® Distal Femur M-O-M

### INSTRUMENTS

#### Content MUTARS® Tibia container 2

**MUTARS® patella drill guide**  
7350-0000



**MUTARS® patella clamp**  
7352-0001



**drill 126 x 3,2 mm**  
4221-0019



**pin inserter 3,2 mm**  
4223-0006



**ic-pin extractor**  
7512-0800



**ic t-handle**  
4223 -0023



**resection check**  
4223 0009



**MUTARS® instrument for locking mechanism**  
7720-1201



**MUTARS® tibial reamer**  
REF 7755-0003



**extractor universal**  
REF 7512-2026



# MUTARS<sup>®</sup> Distal Femur M-O-M

---

## INSTRUMENTS



**setting instrument for locking mechanism**  
7751-1200



**MUTARS<sup>®</sup> positioner for locking mechanism**  
7610-0003



**MUTARS<sup>®</sup> tibial punch**  
7755 0004  
7755-0028 xs



**MUTARS<sup>®</sup> tibia preparator**  
7755-0021



**MUTARS<sup>®</sup> tibia extractor M-O-M**  
7755-0020



**MUTARS<sup>®</sup> reamer for stem preparation**  
7330-1003



**MUTARS<sup>®</sup> assembling forceps**  
7720-1202



## INSTRUMENTS

### Content MUTARS® modular tibia trial container

#### MUTARS® trial for tibial joint M-O-M

7751-0303	extra small
7751-0300	small
7751-0305	standard
7751-0310	large



#### MUTARS® trial tibial spacer

7800-2500	25 mm small
7800-3500	35 mm small
7800-5000	50 mm small



#### MUTARS® trial tibial spacer

7810-0500	5 mm rl lm
7805-0500	5 mm ll rm
7810-1000	10 mm rl lm
7805-1000	10 mm ll rm
7810-1500	15 mm rl lm
7805-1500	15 mm ll rm
7810-2000	20 mm rl lm
7805-2000	20 mm ll rm



#### MUTARS® trial stem for tibial joint

7755-1211	11/120 mm tibial; 11/160mm femoral
7755-1213	13/120 mm tibial; 13/160mm femoral
7755-1215	15/120 mm tibial; 15/160mm femoral
7755-1217	17/120 mm tibial; 17/160mm femoral
7755-1611	11/160 mm tibial; 11/200mm femoral
7755-1613	13/160 mm tibial; 13/200mm femoral
7755-1615	15/160 mm tibial; 15/200mm femoral
7755-1617	17/160 mm tibial; 17/200mm femoral
7755-2011	11/200 mm tibial; 11/240mm femoral
7755-2013	13/200 mm tibial; 13/240mm femoral
7755-2015	15/200 mm tibial; 15/240mm femoral
7755-2017	17/200 mm tibial; 17/240mm femoral



#### MUTARS® trial locking mechanism

7720-1200



#### MUTARS® counter instrument for tibial joint

7755-0027



## INSTRUMENTS

### Content MUTARS® trial components



#### **MUTARS® trial prox. femur**

7710 0205 50 mm

7710 0207 50 mm



#### **MUTARS® trial reducer**

7730 0220 20 mm

7730 0230 30 mm



#### **MUTARS® trial connecting part 100 mm**

7730 0100



#### **MUTARS® trial extension piece**

7772 2504 40 mm

7772 2506 60 mm

7772 2508 80 mm

7772 2510 100 mm



#### **MUTARS® trial extension piece**

77500105 105 mm

77500125 125 mm



#### **MUTARS® trial femoral stem**

7760 0011 11 mm

7760 0013 13 mm

7760 0015 15 mm

7760 0017 17 mm



#### **MUTARS® trial bar screw**

7792 1002 M10x25 mm

7792 1004 M10x45 mm

7792 1006 M10x65 mm

7792 1008 M10x85 mm

7792 1010 M10x105 mm

7792 1012 M10x125 mm

7792 1014 M10x145 mm

7792 1016 M10x165 mm

7792 1018 M10x185 mm

7792 1020 M10x205 mm



## MUTARS® Distal Femur M-O-M

---

### INSTRUMENTS

#### Content MUTARS® rigid drills container

#### MUTARS® rigid drill

4220 4010.1	10 mm
4220 4011.1	11 mm
4220 4012.1	12 mm
4220 4013.1	13 mm
4220 4014.1	14 mm
4220 4015.1	15 mm
4220 4016.1	16 mm
4220 4017.1	17 mm
4220 4018.1	18 mm



**ic adapter outside A/O, inside ic  
canulated**  
7512-3602





## MUTARS<sup>®</sup> Distal Femur M-O-M

---

### INSTRUMENTS

#### **Content MUTARS<sup>®</sup> distal femur M-O-M trial component tray**



#### **MUTARS<sup>®</sup> dist. Femur M-O-M trial implant**

7720-0045	110 mm left
7720-0047	90 mm left
7720-0040	110 mm right
7720-0042	90 mm right

#### **MUTARS<sup>®</sup> trial locking mechanism 7720-1200**



implantcast GmbH  
Lüneburger Schanze 26  
D-21614 Buxtehude  
Germany  
phone: +49 4161 744-0  
fax: +49 4161 744-200  
e-mail: [info@implantcast.de](mailto:info@implantcast.de)  
internet: [www.implantcast.de](http://www.implantcast.de)



Your local distributor:

MUTM2OPE-070813

